

Anti-Cd144/VE Cadherin Monoclonal Antibody For Fc - P33151

Item Number: CM0003987

Introduction

Rabbit monoclonal antibody targeting human CD144/VE Cadherin, conjugated to for direct flow cytometry. Ideal for endothelial cell immunophenotyping and vascular biology research, facilitating adhesion and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD144/VE Cadherin mAb
Remarks/Alias	7B4; CD144
Species	Human
Gene ID (Human)	1003
Gene ID	1003
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 48-593 of human CD144/VE Cadherin(NP_001786.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P33151
Protein Weight	88kDa
Shipping	Ice bag