


Anti-Cd146 Polyclonal Antibody For Wb, Elisa - P43121

Item Number: CM0021329



Introduction

Rabbit polyclonal antibody targeting human CD146 (MUC18) for Western blot and ELISA. Cross-reacts with mouse and rat. Investigate cell adhesion, melanoma metastasis, and endothelial biology.

[Learn More](#)

Parameter	Value
Product Name	CD146 Rabbit pAb
Remarks/Alias	CD146; MUC18; HEMCAM; METCAM; MeICAM
Species	Human
Gene ID (Human)	4162
Gene ID	4162
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 546-645 of human CD146 (NP_006491.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43121
Protein Weight	72kDa
Shipping	Ice bag