

Anti-Cd146 Polyclonal Antibody For Wb, Elisa - Q8R2Y2

Item Number: CM0019677

Introduction

Rabbit polyclonal CD146 antibody validated for WB and ELISA, reactive with human and mouse. Targets cell adhesion molecule MUC18, suitable for melanoma and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	CD146 Rabbit pAb
Remarks/Alias	CD146; CD149; Muc18; s-endo; 1-gicerin; s-gicerin
Species	Mouse
GeneID (Human)	4162
GeneID	84004
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-563 of mouse CD146. (NP_075548.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8R2Y2
Protein Weight	71kDa
Shipping	Ice bag