

Anti-Cd146 Rabbit Monoclonal Antibody For Wb, Ihc-P, Fc, Elisa - Q8R2Y2

Item Number: CM0004699

Introduction

High specificity CD146 Rabbit Monoclonal Antibody validated for WB, IHC-P, FC, ELISA with reactivity across Human, Mouse, Rat. Ideal for cell adhesion and melanoma research.

[Learn More](#)

Parameter	Value
Product Name	CD146 Rabbit mAb
Remarks/Alias	CD146; CD149; Muc18; s-endo; 1-gicerin; s-gicerin
Species	Mouse
GeneID (Human)	4162
GeneID	84004
Immunogen	Recombinant protein corresponding to a sequence within amino acids 24-563 of mouse CD146 (NP_075548.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, FC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8R2Y2
Protein Weight	71kDa
Shipping	Ice bag