


Anti-Cd146/Mcam Monoclonal Antibody For Wb, If/Icc, Elisa - P43121

Item Number: CM0011459



Introduction

High-performance recombinant monoclonal antibody targeting human CD146/MCAM (P43121, 72kDa). Validated for WB, IF/ICC, and ELISA. Plays a role in cell adhesion and tumor metastasis, ideal for cancer and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	CD146/MCAM Rabbit mAb
Remarks/Alias	CD146; MUC18; HEMCAM; METCAM; MelCAM
Species	Human
GeneID	4162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-559 of human CD146/MCAM (NP_006491.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P43121
Protein Weight	72kDa
Shipping	Ice bag