

Anti-Cd146/Mcam Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P43121

Item Number: CM0006916



Introduction

Rabbit monoclonal antibody against human CD146/MCAM (MUC18). Validated for WB, IF/ICC, ELISA. Immunogen: recombinant protein (aa 26-350). Supports vascular biology and melanoma metastasis research. ~72 kDa.

[Learn More](#)

Parameter	Value
Product Name	CD146/MCAM Rabbit mAb
Remarks/Alias	CD146; MUC18; HEMCAM; METCAM; MelCAM
Species	Human
Gene ID (Human)	4162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-350 of human CD146/MCAM (NP_006491.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P43121
Protein Weight	72kDa
Shipping	Ice bag