

Anti-Cd146/Mcam Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P43121

Item Number: CM0007729



Introduction

Highly specific rabbit monoclonal antibody targeting human CD146/MCAM (P43121), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for cell adhesion, melanoma progression, and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	CD146/MCAM Rabbit mAb
Remarks/Alias	CD146; MUC18; HEMCAM; METCAM; MelCAM; CD146/MCAM
Species	Human
Gene ID	4162 (Human)
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 547-646 of human CD146/MCAM (P43121)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43121
Protein Weight	72kDa
Shipping	Ice bag