

Anti-Cd146/Mcam Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc - P43121

Item Number: CM0010982



Introduction

Rabbit monoclonal antibody against human CD146/MCAM (aliases MUC18, MelCAM), validated for WB, IF/ICC, IHC-P, ELISA. Recombinant immunogen, human cross-reactivity; ideal for melanoma progression and endothelial adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	CD146/MCAM Rabbit mAb
Remarks/Alias	CD146; MUC18; HEMCAM; METCAM; MelCAM
Species	Human
GeneID (human)	4162
GeneID	4162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-559 of human CD146/MCAM (NP_006491.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P43121
Protein Weight	72 kDa
Shipping	Ice bag