

Anti-Cd146/Mcam Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P43121

Item Number: CM0016237



Introduction

Rabbit polyclonal antibody against human CD146/MCAM for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Useful for cell adhesion and melanoma metastasis studies.

[Learn More](#)

Parameter	Value
Product Name	CD146/MCAM Rabbit pAb
Remarks/Alias	CD146; MUC18; HEMCAM; METCAM; MelCAM; CD146/MCAM
Species	Human
Gene ID (Human)	4162
Gene ID	4162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-350 of human CD146/MCAM (NP_006491.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43121
Protein Weight	72kDa
Shipping	Ice bag