

Anti-Cd147 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P18572

Item Number: CM0002382



Introduction

Anti-CD147 rabbit monoclonal antibody validated for WB, IHC-P, and ELISA in mouse samples. Targets basigin (EMMPRIN) for retinal development, immune cell adhesion, and monocarboxylate transport studies.

[Learn More](#)

Parameter	Value
Product Name	CD147 Rabbit mAb
Remarks/Alias	HT-7; CD147; EMMPRIN
Species	Mouse
GeneID (Human)	682
GeneID	12215
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 141-302 of mouse CD147 (NP_033898.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	P18572
Protein Weight	42kDa
Shipping	Ice bag