

Anti-Cd147/BSg Rabbit Monoclonal Antibody For Wb And Elisa - P35613

Item Number: CM0006123



Introduction

High-quality rabbit monoclonal antibody against human CD147/BSG (Basigin, EMMPRIN) for Western blot and ELISA. Cross-reacts with mouse; supports research into tumor progression, angiogenesis, and microbial infection.

[Learn More](#)

Parameter	Value
Product Name	CD147/BSG Rabbit mAb
Remarks/Alias	OK; 5F7; TCSF; CD147; EMPRIN; HAB18G; EMMPRIN; CD147/BSG
Species	Human
Gene ID (Human)	682
Gene ID	682
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CD147/BSG (P35613).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P35613
Protein Weight	19kDa/22kDa/29kDa/42kDa
Shipping	Ice bag