

Anti-Cd147/Bsg Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P35613

Item Number: CM0007067



Introduction

Recombinant rabbit monoclonal antibody targeting human CD147/BSG (aa 286-385 epitope) for Western blot, IHC-P, and ELISA. Suitable for cancer biology, retinal development, and host-pathogen interaction studies.

[Learn More](#)

Parameter	Value
Product Name	CD147/BSG Rabbit mAb
Remarks/Alias	OK; 5F7; TCSF; CD147; EMPRIN; HAb18G; EMMPRIN; CD147/BSG
Species	Human
Gene ID (Human)	682
Gene ID	682
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P35613
Protein Weight	42kDa
Shipping	Ice bag