

Anti-Cd147/Bsg Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P35613

Item Number: CM0016940



Introduction

Rabbit polyclonal antibody to human CD147/BSG (Basigin/EMMPRIN) for WB, IHC-P, IP, ELISA. Cross-reacts with mouse and rat. CD147 functions in tumor invasion, angiogenesis, immune cell signaling, and pathogen entry. Suitable for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	CD147/BSG Rabbit pAb
Remarks/Alias	OK; 5F7; TCSF; CD147; EMPRIN; HAB18G; EMMPRIN; CD147/BSG
Species	Human
GeneID (Human)	682
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 286-385 of human CD147/BSG (NP_001719.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35613
Protein Weight	42kDa
Shipping	Ice bag