

Anti-Cd150/Slam Monoclonal Antibody For Fc - Q9Qum4

Item Number: CM0005180

Introduction

Mouse reactive rabbit monoclonal anti-CD150/SLAM antibody, validated for flow cytometry. Targets the 38 kDa signaling lymphocytic activation molecule (SLAMF1) involved in immune cell activation and differentiation. Ideal for immunological research.

[Learn More](#)

Parameter	Value
Product Name	CD150/SLAM Rabbit mAb
Remarks/Alias	Slam; CD150; IPO-3; CDw150; ESTM51; 4933415F16; CD150/SLAM
Species	Mouse
GeneID (human)	6504
GeneID	27218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-242 of mouse CD150/SLAM (NP_038758.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QUM4
Protein Weight	38kDa
Shipping	Ice bag