

Anti-Cd150/Slam Polyclonal Antibody For Wb And Elisa - Q13291

Item Number: CM0022397



Introduction

Rabbit polyclonal antibody targeting human CD150/SLAM (SLAMF1). Suitable for Western blot and ELISA, cross-reacts with human and mouse. SLAM regulates immune cell activation, differentiation, and measles virus entry.

[Learn More](#)

Parameter	Value
Product Name	CD150/SLAM Rabbit pAb
Remarks/Alias	SLAM; CD150; CDw150; CD150/SLAM
Species	Human
GeneID (Human)	6504
GeneID	6504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-237 of human CD150/SLAM (NP_003028.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13291
Protein Weight	37kDa
Shipping	Ice bag