

Anti-Cd151 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P48509

Item Number: CM0024279

Introduction

High-quality rabbit polyclonal antibody against human CD151 (UniProt P48509), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for research on tetraspanin-enriched microdomains, cell adhesion, migration, and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	CD151 Rabbit pAb
Remarks/Alias	EB57; GP27; MER2; RAPH; SFA1; PETA-3; TSPAN24; CD151
Species	Human
GeneID (Human)	977
GeneID	977
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 113-221 of human CD151 (NP_004348.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48509
Protein Weight	28kDa
Shipping	Ice bag