

## Anti-Cd152/Ctla-4 Polyclonal Antibody For Wb, Elisa - P16410

Item Number: CM0013365



### Introduction

Rabbit polyclonal antibody against human CTLA-4 (CD152) validated for WB and ELISA, cross-reacting with mouse and rat. Detects a ~25 kDa band; ideal for T-cell co-inhibition and immuno-oncology research.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD152/CTLA-4 Rabbit pAb                                                                                                  |
| Remarks/Alias    | CD, GSE, GRD4, ALPS5, CD152, CTLA-4, IDDM12, CELIAC3, CD152/CTLA-4                                                       |
| Species          | Human                                                                                                                    |
| GeneID (Human)   | 1493                                                                                                                     |
| GeneID           | 1493                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 49-223 of human CD152/CTLA-4 (NP_005205.2) |
| Source           | Rabbit                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                    |
| Application      | WB, ELISA                                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| SWISS            | P16410                                                                                                                   |
| Protein Weight   | 25kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                  |