

Anti-Cd152/Ctla-4 Rabbit Monoclonal Antibody For Wb And Elisa - P16410

Item Number: CM0007331



Introduction

Rabbit monoclonal antibody targeting human CD152/CTLA-4 (CTLA4), suitable for Western blot and ELISA. Immunogen: recombinant protein. Cross-reacts with human. For research use in T-cell immunoregulation and immune checkpoint studies.

[Learn More](#)

Parameter	Value
Product Name	CD152/CTLA-4 Rabbit mAb
Remarks/Alias	CD; GSE; GRD4; ALPS5; CD152; CTLA-4; IDDM12; CELIAC3; CD152/CTLA-4
Species	Human
GeneID	1493
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-162 of human CD152/CTLA-4 (NP_005205.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P16410
Protein Weight	25kDa
Shipping	Ice bag