

Anti-Cd153 Monoclonal Antibody For Flow Cytometry - P32971

Item Number: CM0004869

Introduction

Human CD153 (CD30L) monoclonal antibody conjugated to , validated for flow cytometry. Ideal for studying T-cell proliferation and CD30/CD30L signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD153 mAb
Remarks/Alias	CD153; CD30L; CD30LG; TNLG3A
Species	Human
GeneID (Human)	944
GeneID	944
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 63-234 of human CD153 (NP_001235.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32971
Protein Weight	26kDa
Shipping	Ice bag