


Anti-Cd155 (Pvr/Necl-5) Rabbit Polyclonal Antibody For Wb And Elisa - P15151

Item Number: CM0018844



Introduction

Anti-CD155 rabbit polyclonal antibody validated for Western blot and ELISA. Recognizes human poliovirus receptor (PVR/NECL-5) involved in NK cell adhesion, tumor invasion, and viral entry. Multiple isoforms ~39-45 kDa.

[Learn More](#)

Parameter	Value
Product Name	CD155 Rabbit pAb
Remarks/Alias	PVR; CD155; HVED; NECL5; Necl-5; PVS; TAGE4; poliovirus receptor
Species	Human
Gene ID (Human)	5817
Gene ID	5817
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15151
Protein Weight	39kDa/40kDa/42kDa/45kDa
Shipping	Ice bag