


Anti-Cd155/Pvr Monoclonal Antibody For Wb, Ihc-P, Elisa - P15151

Item Number: CM0005081



Introduction

CD155/PVR Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human CD155/PVR, a NK cell adhesion molecule and poliovirus receptor, involved in tumor immune evasion. Ideal for cancer and immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD155/PVR Rabbit mAb
Remarks/Alias	PVS; HVED; CD155; NECL5; TAGE4; Necl-5
Species	Human
GeneID (Human)	5817
GeneID	5817
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-343 of human CD155/PVR (NP_006496.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P15151
Protein Weight	45kDa
Shipping	Ice bag