


Anti-Cd155/Pvr Polyclonal Antibody For Wb, Elisa - P15151

Item Number: CM0027506



Introduction

Rabbit polyclonal antibody against human CD155/PVR for Western Blot (WB) and ELISA. Validated cross-reactivity with mouse. Ideal for NK cell adhesion, tumor immunology, and virology research.

[Learn More](#)

Parameter	Value
Product Name	CD155/PVR Rabbit pAb
Remarks/Alias	PVS; HVED; CD155; NECL5; TAGE4; Necl-5; CD155/PVR
Species	Human
Gene ID (Human)	5817
Gene ID	5817
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-345 of human CD155/PVR (NP_006496.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P15151
Protein Weight	45kDa
Shipping	Ice bag