


Anti-Cd157/Bst-1 Monoclonal Antibody For Flow Cytometry - Q10588

Item Number: CM0005532



Introduction

-conjugated rabbit anti-human CD157/BST-1 monoclonal antibody for flow cytometry. Targets CD157/BST-1, a GPI-anchored enzyme involved in cyclic ADP-ribose signaling. Suitable for human sample immunophenotyping and functional studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD157/BST-1 mAb
Remarks/Alias	CD157; cADPR2; BST1
Species	Human
Gene ID	683
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-292 of human CD157/BST-1 (NP_004325.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q10588
Protein Weight	36kDa/29kDa
Shipping	Ice bag