

Anti-Cd16 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08637/O75015

Item Number: CM0007702



Introduction

High-specificity rabbit monoclonal antibody against human CD16 (FCGR3A/FCGR3B), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, rat, and cynomolgus monkey. Ideal for ADCC and NK cell research.

[Learn More](#)

Parameter	Value
Product Name	CD16 Rabbit mAb
Remarks/Alias	CD16; FCG3; CD16A; FCGR3; IGFR3; IMD20; FCR-10; FCRIII; CD16-II; FCGRIII; FCRIIIA; FcGRIIIA
Species	Human
Gene ID (Human)	2214/2215
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-208 of human CD16b (NP_000561.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat, Cynomolgus monkey
SWISS	P08637/O75015
Protein Weight	29 kDa/26 kDa
Shipping	Ice bag