

Anti-Cd161 Rabbit Monoclonal Antibody For Wb, If-P, Fc, Elisa - Q12918

Item Number: CM0006314



Introduction

Rabbit monoclonal antibody against human CD161 (KLRB1). Validated for WB, IF-P, FC, ELISA. Detects 25 kDa protein. Cross-reacts with human, mouse, rat. Ideal for NK cell and immune research.

[Learn More](#)

Parameter	Value
Product Name	CD161 Rabbit mAb
Remarks/Alias	NKR; CD161; CLEC5B; NKR-P1; NKR-P1A; hNKR-P1A; KLRB1
Species	Human
GeneID	3820
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 67-225 of human CD161 (NP_002249.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, FC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12918
Protein Weight	25kDa
Shipping	Ice bag