

Anti-Cd162/Psgl-1 Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - Q14242

Item Number: CM0010780



Introduction

Anti-CD162/PSGL-1 Rabbit Monoclonal Antibody validated for WB, IF/ICC, FC, ELISA in Human. Target is a selectin ligand involved in leukocyte rolling and inflammation.

[Learn More](#)

Parameter	Value
Product Name	CD162/PSGL-1 Rabbit mAb
Remarks/Alias	CLA; CD162; PSGL1; PSGL-1; CD162/PSGL-1
Species	Human
GeneID (Human)	6404
GeneID	6404
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	Q14242
Protein Weight	43 kDa/45 kDa
Shipping	Ice bag