

Anti-Cd162/Psgl-1 Rabbit Monoclonal Antibody For If-P, Ihc-P, Fc, Elisa - Q3Ta56

Item Number: CM0011277

Introduction

Rabbit monoclonal antibody targeting mouse CD162/PSGL-1 (SWISS Q3TA56). Validated for IF-P, IHC-P, FC, and ELISA applications. Cross-reacts with mouse and rat. Ideal for research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	CD162/PSGL-1 Rabbit mAb
Remarks/Alias	CD162; Psgl1; Selp1; Selp; Psgl-1
Species	Mouse
GeneID (Human)	6404
GeneID	20345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-320 of mouse CD162/PSGL-1 (NP_033177.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, FC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q3TA56
Protein Weight	44kDa
Shipping	Ice bag