

Anti-Cd162/PSgl-1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14242

Item Number: CM0007644



Introduction

Rabbit monoclonal antibody against human CD162/PSGL-1, a P-selectin glycoprotein ligand 1. Validated for WB, IHC-P, IF/ICC, and ELISA. Suitable for inflammation and leukocyte adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	CD162/PSGL-1 Rabbit mAb
Remarks/Alias	CLA; CD162; PSGL1; PSGL-1; CD162/PSGL-1
Species	human
GeneID (human)	6404
GeneID	6404
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-295 of human CD162/PSGL-1 (NP_002997.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q14242
Protein Weight	43KDa/44KDa
Shipping	Ice bag