

Anti-Cd162/Psgl-1 Rabbit Polymab® Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - Q14242 / Q3Ta56

Item Number: CM0008081



Introduction

Rabbit monoclonal antibody against human/mouse CD162/PSGL-1, validated for WB, IF/ICC, IF-P, IHC-P, ELISA applications. Detects human, mouse, and rat PSGL-1, a key mediator of leukocyte rolling in inflammation.

[Learn More](#)

Parameter	Value
Product Name	CD162/PSGL-1 Rabbit PolymAb®
Remarks/Alias	CLA; CD162; PSGL1; PSGL-1; CD162/PSGL-1
Species	Human/Mouse
Gene ID (Human)	6404
Gene ID	6404/20345
Immunogen	Recombinant protein of human/mouse CD162/PSGL-1.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14242/Q3TA56
Protein Weight	43kDa/42kDa
Shipping	Ice bag