


Anti-Cd163 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q86Vb7

Item Number: CM0026152



Introduction

Rabbit polyclonal antibody against human CD163, validated for WB, IHC-P, IF/ICC, ELISA. Detects endogenous CD163 in human, mouse, rat samples. Targets the 125 kDa hemoglobin scavenger receptor. Ideal for macrophage and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	CD163 Rabbit pAb
Remarks/Alias	M130; MM130; SCAR11; CD163
Species	Human
GeneID (human)	9332
GeneID	9332
Immunogen	Synthetic peptide/A synthetic peptide corresponding to a sequence within amino acids 750-1050 of human CD163 (NP_004235.4?).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86VB7
Protein Weight	125kDa
Shipping	Ice bag