


Anti-Cd163 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q86Vb7

Item Number: CM0007185



Introduction

Recombinant rabbit monoclonal antibody targeting human CD163 (M130, hemoglobin scavenger receptor) validated for Western blotting (WB), immunofluorescence on paraffin sections (IF-P), and ELISA. Ideal for studying macrophage activation and inflammatory conditions.

[Learn More](#)

Parameter	Specification
Product Name	CD163 Rabbit mAb
Remarks/Alias	M130, MM130, SCAR11, CD163
Species	Human
GeneID (Human)	9332
GeneID	9332
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 861-961 of human CD163 (NP_004235.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	Q86VB7
Protein Weight	125kDa
Shipping	Ice bag