


Anti-Cd163 Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - Q86Vb7 / Q2Vlh6

Item Number: CM0004316



Introduction

High-quality rabbit monoclonal antibody targeting CD163 (M130), validated for WB, IHC-P, IF-P, ELISA in human, mouse, rat, cynomolgus monkey samples. Ideal for macrophage and inflammation research, detecting ~121-125 kDa protein.

[Learn More](#)

Product Name	CD163 Rabbit PolymAb®
Remarks/Alias	M130; MM130; SCAR11; CD163v2; CD163v3
Species	Human/Mouse
GeneID (Human)	9332
GeneID	9332/93671
Immunogen	Recombinant protein of human/mouse CD163.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Cynomolgus monkey
SWISS	Q86VB7/Q2VLH6
Protein Weight	121kDa/125kDa
Shipping	Ice bag