


Anti-Cd163L1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nr16

Item Number: CM0015907



Introduction

Rabbit polyclonal antibody against human CD163L1 (M160), validated for WB and ELISA. Reacts with human and rat. Ideal for studies of scavenger receptor function and macrophage biology.

[Learn More](#)

Parameter	Value
Product Name	CD163L1 Rabbit pAb
Remarks/Alias	WC1; M160; STAR; CD163B; SCAR12; CD163L1
Species	Human
GeneID (Human)	283316
GeneID	283316
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-550 of human CD163L1 (NP_777601.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9NR16
Protein Weight	159kDa
Shipping	Ice bag