


Anti-Cd166/Alcam Monoclonal Antibody () For Flow Cytometry - Q13740

Item Number: CM0012306



Introduction

Recombinant rabbit anti-human CD166/ALCAM monoclonal antibody (clone:) for flow cytometry. Recognizes CD166 (ALCAM), an immunoglobulin superfamily cell adhesion molecule involved in T-cell activation and hematopoiesis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD166/ALCAM mAb
Remarks/Alias	MEMD; CD166
Species	Human
GeneID (Human)	214
GeneID	214
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-526 of human CD166/ALCAM(NP_001618.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13740
Protein Weight	65kDa
Shipping	Ice bag