


Anti-Cd166/Alcam Monoclonal Antibody For Wb, If/Icc, Elisa - Q13740

Item Number: CM0006732



Introduction

Rabbit monoclonal antibody targeting human CD166/ALCAM, validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for studying cell adhesion, T-cell activation, and immunological synapse formation.

[Learn More](#)

Parameter	Value
Product Name	CD166/ALCAM Rabbit mAb
Remarks/Alias	MEMD; CD166; CD166/ALCAM
Species	Human
GeneID (Human)	214
GeneID	214
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-526 of human CD166/ALCAM. (NP_001618.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13740
Protein Weight	65kDa
Shipping	Ice bag