


Anti-Cd166/Alcam Polyclonal Antibody For Wb And Elisa - Q13740

Item Number: CM0021289



Introduction

Rabbit polyclonal antibody targeting human CD166/ALCAM (Q13740). Validated for WB and ELISA, cross-reactive with human and mouse. Ideal for studying cell adhesion, T-cell activation, and hematopoiesis.

[Learn More](#)

Parameter	Value
Product Name	CD166/ALCAM Rabbit pAb
Remarks/Alias	MEMD; CD166; CD166/ALCAM
Species	Human
GeneID (Human)	214
GeneID	214
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-180 of human CD166/CD166/ALCAM (NP_001618.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13740
Protein Weight	65kDa
Shipping	Ice bag