

Anti-Cd168/Rhamm Monoclonal Antibody For Wb, Ihc-P, Elisa - O75330

Item Number: CM0003524

Introduction

Rabbit monoclonal antibody against human CD168/RHAMM (O75330). Ideal for WB, IHC-P, ELISA; cross-reacts with human and mouse. Investigate hyaluronan-mediated cell motility and cancer metastasis.

[Learn More](#)

Parameter	Value
Product Name	CD168/RHAMM Rabbit mAb
Remarks/Alias	CD168; IHABP; RHAMM; CD168/RHAMM
Species	Human
Gene ID (Human)	3161
Gene ID	3161
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75330
Protein Weight	84kDa
Shipping	Ice bag