

Anti-Cd168/Rhamm Polyclonal Antibody For Wb, Elisa - 075330

Item Number: CM0028334

Introduction

Rabbit polyclonal antibody against human CD168/RHAMM, validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for studying hyaluronan-mediated cell motility and metastasis.

[Learn More](#)

Parameter	Value
Product Name	CD168/RHAMM Rabbit pAb
Remarks/Alias	CD168; IHABP; RHAMM; CD168/RHAMM
Species	Human
GeneID (Human)	3161
GeneID	3161
Immunogen	Recombinant protein corresponding to amino acids 341-440 of human CD168/RHAMM (NP_001136028.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75330
Protein Weight	84kDa
Shipping	Ice bag