

Anti-Cd16B Rabbit Monoclonal Antibody For Wb, Fc, Elisa - O75015

Item Number: CM0007904



Introduction

Rabbit monoclonal antibody against human CD16b (FcγRIIIB). Validated for WB, FC, ELISA; detects a 26kDa band. Ideal for studying Fc receptor-mediated immune complex trapping and neutrophil biology.

[Learn More](#)

Parameter	Value
Product Name	CD16b Rabbit mAb
Remarks/Alias	CD16; FCG3; CD16A; CD16b; FCGR3; CD16-I; FCGR3A; FCR-10; FCRIII; FCRIIb
Species	Human
Gene ID (Human)	2215
Gene ID	2215
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 20-208 of human CD16b (NP_000561.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	O75015
Protein Weight	26kDa
Shipping	Ice bag