


Anti-Cd170/Siglec-5 Monoclonal Antibody For Flow Cytometry - O15389

Item Number: CM0005651



Introduction

Rabbit monoclonal antibody against human CD170/Siglec-5, validated for flow cytometry (FC). Suitable for studying myeloid cell adhesion and sialic acid binding in immunological research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD170/Siglec-5 mAb
Remarks/Alias	CD170; OBBP2; CD33L2; OB-BP2; SIGLEC-5
Species	Human
GeneID (Human)	8778
GeneID	8778
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O15389
Protein Weight	61kDa
Shipping	Ice bag