

Anti-Cd171/L1Cam Monoclonal Antibody For Wb, Ihc-P, Elisa - P32004

Item Number: CM0007117



Introduction

Rabbit monoclonal antibody against human CD171/L1CAM, validated for Western blot, IHC (paraffin), and ELISA. Suitable for studying neuronal development, migration, and synaptic plasticity in brain research.

[Learn More](#)

Parameter	Value
Product Name	CD171/L1CAM Rabbit mAb
Remarks/Alias	S10; HSAS; HYCX; MASA; MIC5; SPG1; CAML1; CD171; HSAS1; N-CAML1; NCAM-L1; N-CAM-L1; CD171/L1CAM
Species	Human
GeneID (Human)	3897
GeneID	3897
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P32004
Protein Weight	140kDa
Shipping	Ice bag