

Anti-Cd171/L1Cam Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P32004

Item Number: CM0027826



Introduction

High-quality rabbit polyclonal antibody against CD171/L1CAM, validated for Western blot, IF/ICC, and ELISA. Recognizes human, mouse, and rat L1CAM. Ideal for neuroscience research on neuronal adhesion, migration, and synaptic plasticity.

[Learn More](#)

Parameter	Value
Product Name	CD171/L1CAM Rabbit pAb
Remarks/Alias	S10; HSAS; HYCX; MASA; MIC5; SPG1; CAML1; CD171; HSAS1; N-CAML1; NCAM-L1; N-CAM-L1; CD171/L1CAM
Species	Human
Gene ID (Human)	3897
Gene ID	3897
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1143-1257 of human CD171/L1CAM (NP_000416.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32004
Protein Weight	140kDa
Shipping	Ice bag