


Anti-Cd172A/Sirpα Polyclonal Antibody For Wb, If/Icc, Elisa - P78324

Item Number: CM0026579



Introduction

CamelBio offers anti-CD172a/SIRPα rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat SIRPα (55 kDa). Ideal for studying phagocytosis regulation, immune signaling, and synaptic function research.

[Learn More](#)

Parameter	Value
Product Name	CD172a/SIRPα Rabbit pAb
Remarks/Alias	BIT; MFR; P84; SIRP; MYD-1; SHPS1; CD172A; PTPNS1; CD172a/SIRPα
Species	Human
GeneID (human)	140885
GeneID	140885
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-375 of human CD172a/SIRPα (NP_542970.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78324
Protein Weight	55 kDa
Shipping	Ice bag