

Anti-Cd177 Rabbit Monoclonal Antibody For Fc, Elisa - Q8N6Q3

Item Number: CM0012277

Introduction

Recombinant rabbit monoclonal antibody to human CD177 antigen (NB1/PRV-1). Validated for flow cytometry and ELISA, suitable for neutrophil activation, adhesion, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	CD177 Rabbit mAb
Remarks/Alias	CD177; HNA-2a; HNA2A; NB1; NB1 GP; PRV-1; PRV1; NB1GP; CD177 antigen
Species	Human
GeneID (human)	57126
GeneID	57126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-407 of human CD177 (NP_065139.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Human
SWISS	Q8N6Q3
Protein Weight	15kDa/27kDa/46kDa
Shipping	Ice bag