

Anti-Cd177 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q8N6Q3

Item Number: CM0007132

Introduction

CD177 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Detects human CD177 antigen (NB1/PRV-1/HNA-2a), a neutrophil-specific GPI-anchored protein involved in adhesion and migration. Ideal for inflammation, hematology, and autoimmune disease research.

[Learn More](#)

Parameter	Value
Product Name	CD177 Rabbit mAb
Remarks/Alias	NB1; PRV1; HNA2A; PRV-1; HNA-2a; NB1 GP; CD177
Species	Human
Gene ID (Human)	57126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-407 of human CD177(NP_065139.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q8N6Q3
Protein Weight	51kDa
Shipping	Ice bag