

Anti-Cd18 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P05107

Item Number: CM0008572



Introduction

Rabbit monoclonal anti-human CD18 (ITGB2) antibody validated for WB, IF-P, IHC-P, ELISA. Recognizes human, mouse, rat integrin beta-2 (85 kDa), ideal for leukocyte adhesion and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	CD18 Rabbit mAb
Remarks/Alias	LAD; CD18; MF17; MFI7; LCAMB; LFA-1; MAC-1
Species	Human
GeneID	3689
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD18 (NP_000202.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05107
Protein Weight	85kDa
Shipping	Ice bag