

Anti-Cd19 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P15391

Item Number: CM0004010



Introduction

Rabbit monoclonal antibody targeting human CD19. Suitable for WB, IF-P, IHC-P, and ELISA. Crucial for B-cell receptor signaling and immunodeficiency research.

[Learn More](#)

Parameter	Value
Product Name	CD19 Rabbit mAb
Remarks/Alias	B4; CVID3; CD19
Species	Human
GeneID (Human)	930
GeneID	930
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 456-556 of human CD19. (NP_001761.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P15391
Protein Weight	61kDa
Shipping	Ice bag