

Anti-Cd19 Polyclonal Antibody For Wb, If/Icc, Elisa - P15391

Item Number: CM0024854



Introduction

CD19 rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Recognizes human, mouse, and rat CD19. Ideal for B-cell signaling studies and lymphocytic leukemia research. Your one-stop IVD raw-material partner.

[Learn More](#)

Parameter	Value
Product Name	CD19 Rabbit pAb
Remarks/Alias	B4; CVID3; CD19
Species	Human
Gene ID (Human)	930
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 457-556 of human CD19 (NP_001761.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15391
Protein Weight	61kDa
Shipping	Ice bag