

Anti-Cd19 Rabbit Monoclonal Antibody For Wb, Elisa - P15391

Item Number: CM0006763



Introduction

Rabbit monoclonal antibody against human CD19 protein. Suitable for western blot (WB) and ELISA applications. Reacts with human, mouse, and rat CD19. Ideal for B-cell immunology research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	CD19 Rabbit mAb
Remarks/Alias	B4; CVID3; CD19
Species	Human
GeneID (Human)	930
Immunogen	Recombinant protein of human CD19
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15391
Protein Weight	61kDa
Shipping	Ice bag