


Anti-Cd1B Rabbit Monoclonal Antibody For If/Icc, Fc, Elisa - P29016

Item Number: CM0008012



Introduction

Rabbit monoclonal antibody targeting human CD1B. Validated for IF/ICC, FC, and ELISA. Detects CD1B isoforms (37kDa/30kDa) involved in lipid antigen presentation and T-cell immunity. Ideal for immunology studies.

[Learn More](#)

Parameter	Specification
Product Name	CD1B Rabbit mAb
Remarks/Alias	R1; CD1; CD1A
Species	Human
GeneID (human)	910
GeneID	910
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-303 of human CD1B (NP_001755.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P29016
Protein Weight	37kDa/30kDa
Shipping	Ice bag